

RITRAMA

Ri-Mark Transparent - Colour definition



Light Yellow
Transparent 706
LAB: 80.38/19.7/91.23
RGB: 255/183/0
CMYK: 0/31/100



Yellow
Transparent 710
LAB: 75.47/32.86/80.92
RGB: 255/159/0
CMYK: 0/45/100



Orange
Transparent 720
LAB: 65.79/54.65/57.75
RGB: 255/113/57
CMYK: 0/70/82



Light Red
Transparent 729 Light Red
LAB: 49.13/62.73/41.97
RGB: 214/51/50
CMYK: 10/95/90



Red
Transparent 730
LAB: 45.76/59.4/35.5
RGB: 200/50/53
CMYK: 15/95/86



Magenta
Transparent 789
LAB: 48.93/63.45/-0.99
RGB: 206/54/121
CMYK: 16/93/24



Brown
Transparent 745
LAB: 32.59/18.98/11.16
RGB: 109/66/61
CMYK: 42/71/66



Blue
Transparent 750
LAB: 45.85/-24.13/-53.13
RGB: 0/123/198
CMYK: 100/41/3



Dark Blue
Transparent 752
LAB: 38.05/0.95/-56.91
RGB: 0/93/183
CMYK: 97/67/0



Blue-Violet
Transparent 740
LAB: 32.44/17.64/-58.18
RGB: 38/69/169
CMYK: 93/83/0



Apple Green
Transparent 778
LAB: 69.62/-33.89/57.7
RGB: 131/187/52
CMYK: 54/5/100



Light Green
Transparent 767
LAB: 52.13/-63.25/23.06
RGB: 0/148/81
CMYK: 99/13/95



Green
Transparent 770
LAB: 49.82/-48.86/35.06
RGB: 0/138/53
CMYK: 88/20/100



Grey
Transparent 704
LAB: 63.24/-1.13/7.33
RGB: 155/153/140
CMYK: 42/34/44